

## Oxyhydrogen Ampoule Bottle Sealing Machine



### Okay Energy Equipment Co.,Ltd

http://www.okayenergy.com/





### (1) Introduction of Oxy-hydrogen sealing machine

Traditional Ampoule bottle sealing adopt liquefied gas as fuel, use air as combustion-supporting, so its flame have CO2, CO and nitride which may have bad influence to the medicine in the ampoule bottle, and it is also not safe for working environment.

Okay Energy Oxy-hydrogen Generator have huge Success on Ampoule bottle sealing .Oxy-hydrogen generator mainly electrolyzes water to hydrogen and oxygen, when burn the gas, the oxy-hydrogen flame temperature is up to 2800°C, the flame can seal ampoule bottle rapidly. And the effect as well as the speed is higher and faster than adopting by liquefied gas.

# (2) Comparison of oxy-hydrogen flame sealing with traditional liquefied gas sealing method.

Advantages by Oxy-hydrogen flame sealing on ampoule bottle:

1. Safety: Oxy-hydrogen technology only need water and a little electricity, no need extra cylinder to store gas, it is plug and play, user have no worry about gas leakage or explosion.

2. Convenience: Oxy-hydrogen sealing only need water and a little electricity. One liter water can be used for 7-8 hours. Traditional liquefied gas sealing need cylinder store gas which need be moved, inconvenient.

3. Clean: When Oxy-hydrogen gas burns, it can bring very little invisible water va por, which have no harm to medicine in ampoule bottle. However, when liquefied gas b urns, there is Co2, Co and NOx which may influence the medicine quality in the bottle.

4. High efficiency: The speed to seal ampoule bottle by Oxy-hydrogen is faster tha n that by liquefied gas.

#### (3) Composition of Oxy-hydrogen flame sealing

Oxy-hydrogen generator +torch + adjustable frame supporting + tweezers.





Tel:86-731-84169262 Mail:sherryhe@okayenergy.cn Web:www.okayenergy.com




## (4) How to operate the oxy-hydrogen generator for sealing on ampoule bottle?

1. Install the oxy-hydrogen generator by following manual, then connect the gas pipe, and adjust the position of the torch.

2. Turn on the Oxy-hydrogen generator, the machine start working to produce gas, when the pressure of the piezometer rises up to 0.1Mpa, open the valve of the torch for gas out.



3. Adopt lighter to fire the oxy-hydrogen gas for flame, user can adjust the valve to control the flame.(Attention: when fire the gas, don't let any parts of your body in the front of the flame nozzle in case burnt ;when use lighter , user need fire from the side of the flame torch in case burnt your hand.







Tel:86-731-84169262 Mail:sherryhe@okayenergy.cn Web:www.okayenergy.com

4. There is a tweezers with the machine, please use this professional tweezers to clamp the ampoule bottle tightly in case the bottle falling off, at the same time the bottle should be rotated by hands.



5. Put the bottle which clamped by the tweezers to the oxy-hydrogen gas flame for heating .Rotate the bottle while heating .When the bottle is heat ,pull up the bottle head by tweezers, then bottle is sealed. Please see the video by the link <a href="http://www.okayenergy.com/video/product-video/65.html">http://www.okayenergy.com/video/product-video/65.html</a>

http://www.okayenergy.com/video/product-video/ampoule-sealing-manufacturing-video.html

